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Intended use

Dental splinting material long-term

Extended purpose

Tough-elastic, highly transparent thermoforming foil for the production of occlusal splints such as orthodontic devices or aligner and retainer splints.

Indications for use

Suited for interim splints in case of myoarthropathy (muscular and mandibular discomfort) and accompanying dysfunctions resulting from it, like pressing, bruxism, etc. The splint provides for a very good bite lock and does not cause any tension. Furthermore the production is very easy and there is no need for anchorages, the splint fits very well and secure due to the snap-in effect. Biolon bonds with acrylate.

Contraindication

In case of allergies to any of the ingredients do not use Biolon. Not desired side effects of this medical device are not to be expected in case of proper processing and application. However, immune reactions (e. g. allergies) or irritations principally cannot be excluded. Should you get to know non-desired side effects – also in doubtful cases – we appreciate in any case to be informed with a detailed description of the attendant circumstances and symptoms. We will follow each hint.

Target patient group

Patients who need dental care.

Intended users

Trained staff in dental field

Processing

All Biolon materials are delivered with an insulating foil. This extremely thin foil (0.07 mm) is thermoformed together with the Biolon on the model side and guarantees a decisively improved transparency of the splint. Thus an additional insulation of the model surface is no longer necessary.

If processed in a Vacfomat, the insulation foil must be removed before thermoforming.

The adjustment of the heating and cooling times can either be taken from the survey of thermoforming foils, or can be

programmed directly into the unit with the help of the barcode on the Biolon packaging.

If a Drufoamat scan or a Drufosmart scan is used, please take note of the special information included in this package. We recommend Cenit L (REF D3815) for the necessary blocking out works that are done on the model. The splints produced from Biolon can be built up with dental acrylic resins.

Cleaning

For cleaning, the use of MaxiFresh™ cleaning tablets is recommended.

Please be aware of the MaxiFresh™ instruction for use.

Note

While working out on the material, it is recommended to use an exhaust system and protective glass or goggles as well as a face mask.

Serious incidents

All serious incidents occurring in relation to the product shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is settled.

Storage conditions

At room temperature, i. e. 18–28 °C.

Protect against humidity, always close foil bag when packs are in use.

Lot number / Expiry date

The lot number and the expiry date are indicates on each Biolon packaging. In case of claims, please always indicate the lot number of the product. Do not use the product after the expiry date.

Disposal

Disposal of contents / container in accordance with local / regional / national / international regulations.

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