

Bio Bleach hard



Technical Data Sheet

Product description: Dental thermoforming blank for the production of stable bleaching/medication splints ending at the gingival margin.

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	Chemical Denomination
1	polyethylene terephthalate (PET)

Physical characteristics:

Colour: transparent

Thickness: 0.5 mm

Diameter: Ø 120mm (special diameters on demand)

Density: 1.27g/cm³ DIN 53479

Elongation at break: > 200% DIN 53455

Rockwell hardness: R104- 117 ASTM D785

Flexural strength: 50N/mm² DIN 53455

Absorption of humidity: 0.2% 24 h at 23°C / 3mm

Temperature range: - 40°C- + 65°C

Dreve Dentamid GmbH

Max-Planck-Str. 31

59423 Unna / Germany

Tel.: 02303 / 8807-0

Fax: 02303 / 8807-55

info@dreve.de

www.dreve.com

These data result from measurements which are executed from time to time
in view of our quality management system.
This document is valid without any signature.

Rev. 00/03.2011
Seite 1 von 1
41-4294/LW